



ULTRASONIC SOLUTION: GENERAL PURPOSE

Advanced formulation combines nonionic surfactants, powerful chelators, and corrosion inhibitors designed to enhance the ultrasonic cleaning action and maintain delicate instruments.

Used with an ultrasonic cleaning unit to remove biomatter, debris and other surface soils from instruments and appliances.

May be used directly in the ultrasonic tank or as a sonic conductor when using Optimize's Tartar & Stain Remover.



PRODUCT	SIZE	ITEM NUMBER
Liquid Concentrate	128 fl oz (3.78 L)	00122
Powder Concentrate	4.23 oz (120 g)	00158

Liquid Concentrate makes 10 gallons (37.85 L) of solution.
Powder Concentrate makes 1 gallon (3.78 L) of solution.

keystone
industries

dental.keystoneindustries.com

US: 480 South Democrat Road Gibbstown NJ, 08027
EU: Keystone Industries GmbH Stockholzstr. 11 78224 Singen, Germany

US: 1.800.333.3131
EU: +49 77 31 91 21 01