



Aegis® Pit & Fissure Sealant with ACP

Aegis® Pit and Fissure Sealant is a light-cured sealant that contains the “smart material” Amorphous Calcium Phosphate (ACP).

About Aegis Pit & Fissure Sealant

The ACP filler is resilient and flexible, creating a stronger and longer lasting sealant. Aegis P&F can be used with any light-curing device and offers the added convenience of dispensing through syringe tips. Aegis' controlled flowability keeps the sealant on the tooth while completely filling the pits and fissures and can be placed on the buccal. When properly maintained and placed, Aegis forms a chemical and thermal barrier that can help keep patients free from carious lesions on the occlusal surfaces. The new resin formula allows for better flow and handling. Sealant cures to an opaque film on enamel.



Features & Benefits

- 2nd generation gear design
- 38% ACP filled, promoting the remineralization of tooth structure
- Resilient & flexible, resulting in increased durability & retention
- Non-reliant on patient compliance
- Opaque in color
- Light-cured

PRODUCT

ITEM NUMBER

Standard Kit

0921348

Standard Kit contains:
(enough material for approximately 30 patients)

- 3 x 1.2 ml Syringes Aegis P&F Sealant
- 1 x 3 ml syringe Acid Etch 35%
- 20 x 25 gauge Flowthru® Microbrush®-style tips
- 20 Grey Dispensing Tips, 22 gauge
- 20 Placement Brushes

Refill Kit contains:

- 4 x 1.2ml Aegis P&F Sealant
- 20 Grey Dispensing Tips, 22 gauge
- 20 Placement Brushes

Microbrush® and Flowthru® are registered trademarks of Microbrush Corporation.

keystone[®]
industries

dental.keystoneindustries.com