



Diamond D[®] Denture Acrylic

Advanced natural-looking dentures require artistry as well as technical expertise.

About Diamond D[®]

Since its introduction in 2003, Diamond D[®] has become one of the world's leading high-impact denture acrylic. Laboratories and dentists alike have come to rely on Diamond D's superior strength and natural fibered colors to deliver a premium, fracture-resistant denture product to patients.

Features and Benefits


- Available in industry standard colors
- High-impact strength
- Exceptional flexural strength
- Consistent & forgiving
- Available in heat cure & self cure*

*Self Cure not for sale in the EU

IMAGE	COLOR	POUNDS	HEAT CURE
	ORIGINAL	1 LB	1013061
		5 LBS	1013065
		25 LBS	1013069
	LIGHT	1 LB	1013062
		5 LBS	1013066
		25 LBS	1013070
	LIGHT REDDISH PINK	1 LB	1013067
		5 LBS	1013067
		25 LBS	1013071
	CHROMA ESSENCE	1 LB	1013064
		5 LBS	1013068
		25 LBS	1013072
	DARK VEINED	1 LB	1013112
		5 LBS	1013113
		25 LBS	1013114

Dental Acrylic Izod Impact Strength

Impact strength - a measure of toughness. The stress to fracture a notched specimen with a single blow. Expressed in foot-pounds of energy absorbed. Designed as "Charpy" or "Izod" impact strength depending on the testing machine used.



Source: Notched Izod Impact Strength Determinations per ASTM D256. Submitted in ASTM D638, Type 1 specimen form. Performed at Bodycote Polymer-Broutman Laboratory, a UL and IAPMO approved independent testing laboratory, A2LA and IAS accredited. BBI Test Reference #: 30002570 and 30000666.



dental.keystoneindustries.com